

FYS2150

Lab Report: Drag

Nicholas Karlsen

April 25, 2018

Abstract

A study on the flow of an assortment of spheres in a fluid and the use of image processing to determine the terminal velocity.

1 Introduction

2 Theory

3 Experimental Procedure

3.1 Video capture

This report contains the description and analysis of data collected in the lab 21.03.2018 concerning the flow of several spherical objects in a large range of different sizes and densities. The balls were immersed in fluid, dropped and filmed. Post-lab, the raw footage was then processed using a Python script in order to quantify the motion of the spheres. This



Figure 1: Signal used to determine the FPS of the camera

In order to capture the motion of the balls, a USB video camera was connected to a

computer running uEye cockpit [1], which we used to change the settings of the camera and make recordings.

First, the error of the stated FPS of the camera had to be determined. This was done by connecting a series of circular LEDs (see Fig. 1) to a signal generator. The light emitted would "circle" at a rate which could be changed using the signal generator. By adjusting the signal such that the emitted light would seem stationary when observed through the video feed in uEye cockpit, the FPS of the camera determined to be




Figure 3: Most of the balls used in the experiment, excluding the ones labeled small 1 and 2.

4 Results

Type	Mass [g]	Diameter [mm]	FPS	T [°C]	Filename
Metal	502.76	48.98	100	22.7	A1.avi
Metal	28.13	19.02	100	22.8	A2.avi
Metal	6.99	11.97	100	22.6	A3.avi
Metal	2.08	7.99	100	22.6	A4.avi
Metal	0.68	5.48	100	22.5	A5.avi
Metal	0.10	2.98	100	22.6	A6.avi
Plastic	488.41	99.4	100	22.5	B1.avi
Plastic	61.56	50.02	100	22.5	B2.avi
Plastic	7.12	23.89	100	22.6	B3.avi
Plastic	0.87	12.06	100	22.6	B4.avi
White	29.74	25.24	100	22.6	C1.avi
BigB lack	31.42	21.08	100	22.5	C4.avi
Small Black	5.67	16.45	100	22.5	C5.avi
Big Green	31.60	21.86	100	22.4	C3.avi
Small Green	5.60	16.38	100	22.3	C6.avi
Big Red	18.44	24.01	100	22.3	C2.avi
Glass	0.27	5.81	100	22.3	C7.avi
Small 2	12.0E-3	1.59	100	23.7	D2.avi
Small 1	4.1E-3	1.0	100	23.7	D1.avi

Table 1: Spheres

5 Discussion

6 Conclusion

References

- [1] DS Imaging Development Systems GmbH. ueye cockpit.

*

A Code

Following

scripts/lesVideo_conv.py

```
1 #!/usr/bin/env python
2 # -*- coding: utf-8 -*-
3 '''
4 Reads video file and converts to binary image
5 resulting in easy data analysis.
6 author: Nicholas Karlsen
7
8 Note: skvideo is not included in anaconda python,
9 install by 'pip install sk-video' in terminal.
10 '''
11
12 import numpy as np
13 import skvideo.io
14 import inspect
15 import os
16 import matplotlib.pyplot as plt
17 from skimage.measure import regionprops
18 from matplotlib.image import imread
19 from skimage import util
20 import FYS2150lib as fys # Used for linfit
21 # import skimage.color
22 # from PIL import Image
23 # import skimage.morphology as morph
24 # from skimage import filters
25
26
27 def rgb2gray(rgb):
28     '''
29     Converts shape=(N,M,rgb) array to (N, M) grayscale array see wiki page
30     '''
31     return np.dot(rgb[..., :3], [0.299, 0.587, 0.114]).astype(int)
32
33
34 def gray2binary(gray, limBW=128):
35     """Converts grayscale image to binary grayscale of 0 OR 255
36     image must be array of shape=(N, M)
37     gray: (N, M) array
38     limBW:
39     """
40     bw = np.asarray(gray).copy()
41     bw[bw < limBW] = 0 # Black
42     bw[bw >= limBW] = 255 # White
43     return bw
44
45 def genFilter(image):
```

```

46     """
47     Generates an array to filter out
48     static background based on first frame
49
50     NOT YET IMPLEMENTED
51     """
52     gsImage = rgb2gray(image)
53     bwImage = gray2binary(gsImage)
54     bwImage = bwImage / 255.0
55     return bwImage.astype(int)
56
57 def trackCircle(filename="litenmetallkule.avi", path="current",
58                 hMin=0, hMax=-1, wMin=0, wMax=-1):
59     """
60     Takes video file as input, filters out static background based on
61     first frame and finds the CM of circle in every frame. Requires
62     circle to be only object in frame (after filtering), so requires static
63     background. If not, try adjust hMin, hMax, wMin, wMax to crop out
64     moving
65     background.
66     filename: filename of video
67     path: FULL path of file, eg '/home/nick/Videos/fys2150drag'.
68           if left as default, it will assume same path as script.
69     hMin, hMax, wMin, wMax: used for cropping the image.
70     """
71
72     # Fetching current dir path
73     folderPath = os.path.dirname(
74         os.path.abspath(
75             inspect.getfile(
76                 inspect.currentframe()))))
77
78     # If path is specified, use that instead.
79     if path != "current":
80         folderPath = path
81         "if path is specified"
82
83     fullFilename = folderPath + "/" + filename
84
85     print "Reading video..."
86
87     video = skvideo.io.vread(fullFilename)
88     totalFrames = len(video)
89
90     print "Number of frames:", len(video)
91
92     frameStart = 0
93     frameStop = totalFrames
94
95     "Creates array to store x, y vals of CM"
96     cmPos = np.zeros([frameStop - frameStart, 2])
97
98     validFrames = [] # Keeps track of usable frames

```

```

98
99 def detectCirc(image):
100     """
101     Inverts color of image and detects center of circle shape.
102     Assumes circle is the ONLY object in image, so noise
103     needs to be filtered out
104     """
105     #staticBg = genFilter(video[0])
106
107     invFrame = image
108     bwFrame = gray2binary(
109         rgb2gray(
110             util.invert(invFrame))) [hMin:hMax, wMin:wMax]
111     bwFrame = bwFrame # * staticBg
112     # Detects shapes in image
113     props = regionprops(label_image=bwFrame.astype(int))
114
115     return props, invFrame, bwFrame
116
117 for frame in xrange(totalFrames):
118     """
119     Need to invert image for regionprops to work, only finds white obj
120     on black background, not black on white.
121     """
122     # convert to binary grayscale to filter out noise
123     props = detectCirc(video[frame])[0]
124
125     # Bad way of checking if the ball is in frame
126     if len(props) == 0:
127         cmPos[frame] = "nan"
128     else:
129         cmPos[frame] = props[0].centroid # Detects centroids
130         validFrames.append(frame) # Keeps track of frames with ball
131
132         # Print info to terminal while processing
133         print "frame", frame, \
134             "_", "Center of mass:", \
135             "x=%i, y=%i" % (cmPos[frame][1], cmPos[frame][0])
136
137 def plot_im(frame=int(totalFrames / 2.0)):
138     "plot frame + CM, used to check functionality"
139     im = video[frame]
140     props, invFrame, bwFrame = detectCirc(im)
141     cmPos[frame] = props[0].centroid
142     plt.subplot(311)
143     plt.imshow(invFrame)
144     plt.title("Raw image, frame:%i" % frame)
145     plt.plot(cmPos[frame, 1], cmPos[frame, 0] + hMin,
146             "ro", label="Center of mass")
147     plt.legend()
148     plt.subplot(312)
149     plt.imshow(bwFrame, cmap=plt.get_cmap('gray'))
150     plt.plot(cmPos[frame, 1], cmPos[frame, 0],

```

```

151         "ro", label="Center of mass")
152     plt.title("Processed image, frame:%i" % frame)
153     plt.legend()
154     plt.show()
155 plot_im()
156
157     return cmPos.astype(int), validFrames
158
159
160 def testFunc():
161     """
162     Testing that method of finding C.M works properly
163     """
164     import skimage.color
165     #img = imread("bilde5.png")
166     img = imread("frame_inv2.png")
167     bwImg = skimage.color.rgb2gray(img)
168     plt.subplot(211)
169     plt.imshow(bwImg, cmap=plt.get_cmap('gray'))
170
171     props = regionprops(label_image=bwImg.astype(int))
172     cm = props[0].centroid
173
174     plt.subplot(212)
175     plt.imshow(img)
176     plt.plot(cm[1], cm[0], "ro",
177             label="Center of mass = (%i, %i)" % (cm[1], cm[0]))
178     plt.legend()
179     plt.show()
180
181
182
183
184 if __name__ == "__main__":
185
186     def readlabdat(filename):
187         """
188         Used to read the file which stores the parameters of the
189         sphere
190         """
191         vids = []; mass = []; radius = []; temp = []
192
193         file = open(filename, "r")
194         for line in file:
195             cols = line.split()
196             mass.append(cols[1])
197             radius.append(cols[2])
198             temp.append(cols[-2])
199             vids.append(cols[-1])
200         file.close()
201
202         return mass, radius, temp, vids
203     mass, radius, temp, vids = readlabdat("data/labdata.dat")

```

```

204
205 folderPath = "/home/nick/Videos/fys2150drag"
206
207 rows = range(1, 19)
208
209 outfile = open("data/results.dat", "w")
210
211 for row in rows:
212     cm, validFrames = trackCircle(filename=str(vids[row]),
213                                   path=folderPath,
214                                   hMin=67, hMax=216)
215
216     if len(cm[:, 1]) != len(validFrames):
217         print "Tracking interrupted in some frames,"
218         print "Only returning uninterrupted frames."
219         x = []
220         y = []
221         for validFrame in validFrames:
222             x.append(cm[validFrame, 1])
223             y.append(cm[validFrame, 0])
224         x = np.array(x).astype(int)
225         y = np.array(y).astype(int)
226     else:
227         x = cm[validFrames[0]:validFrames[-1], 1]
228         y = cm[validFrames[0]:validFrames[-1], 0]
229
230     x = np.array(x)
231     y = np.array(y)
232     validFrames = np.array(validFrames)
233
234     print "Find start/stop of terminal velocity (straight, steep line)
to perform linfit:"
235
236     plt.plot(validFrames, x, "o")
237     plt.xlabel("Frame")
238     plt.ylabel("x-position of center of mass [px]")
239     plt.title("Use to determine start/stop frame of linfit")
240     plt.show()
241
242     start = int(input("Start index:"))
243     stop = int(input("Stop index:"))
244
245     m, c, dm, dc = fys.linfit(validFrames[start:stop], x[start:stop])
246
247     plt.subplot(211)
248     plt.plot(validFrames, x, ".", label="Position of CM")
249     plt.plot(validFrames[start:stop],
250             validFrames[start:stop] * m + c,
251             label="linear fit, y=mx+c")
252     plt.text(0, 1000, "m = %i [px/frame]\n dm = %i [px/frame]" % (m, dm)
253 ))
254     plt.xlabel("Frame")
255     plt.ylabel("x-pos [px]")

```



```

255     plt.legend()
256     plt.subplot(212)
257     plt.plot(validFrames, y, ".", label="Position of CM")
258     plt.xlabel("Frame")
259     plt.ylabel("y-pos [px]")
260     plt.show()
261
262     outfile.write(vids[row] + " & " + "%i"%(m) + " & " + "%i"%dm + "\\n")
263
264     outfile.close()

```

scripts/data/labdata.dat

	Type	Mass [g]	Diameter [mm]	FPS	T[C]	Filename
1						
2	Metal	502.76	48.98	100	22.7	A1.avi
3	Metal	28.13	19.02	100	22.8	A2.avi
4	Metal	6.99	11.97	100	22.6	A3.avi
5	Metal	2.08	7.99	100	22.6	A4.avi
6	Metal	0.68	5.48	100	22.5	A5.avi
7	Metal	0.10	2.98	100	22.6	A6.avi
8	Plastic	488.41	99.4	100	22.5	B1.avi
9	Plastic	61.56	50.02	100	22.5	B2.avi
10	Plastic	7.12	23.89	100	22.6	B3.avi
11	Plastic	0.87	12.06	100	22.6	B4.avi
12	White	29.74	25.24	100	22.6	C1.avi
13	BigBlack	31.42	21.08	100	22.5	C4.avi
14	SmallBlack	5.67	16.45	100	22.5	C5.avi
15	BigGreen	31.60	21.86	100	22.4	C3.avi
16	SmallGreen	5.60	16.38	100	22.3	C6.avi
17	BigRed	18.44	24.01	100	22.3	C2.avi
18	Glass	0.27	5.81	100	22.3	C7.avi
19	Small 1	4.1e-3	1.0	100	23.7	D1.avi
20	Small 2	12.0e-3	1.59	100	23.7	D2.avi